

## Radiant Heater RH600 Technical Service Bulletin # 6

**Re:** List of useable Logic PCBs

**Date:** October 16, 2001

**Reference Doc:** N/A

**Reason:** For your information the list of useable and not useable Logic PCBs. Please note that all devices RH600 have been updated during the actions according to TSB RH600 # 5 with updated and useable Logic PCBs.

**Solution:** Table of Logic PCBs and status:

P/N Logic PCB	Status
8306101 / SE8306101	Not useable
2M22408 / SE2M22408	Useable
2M22409 / SE2M22409	Useable
2M22383 / SE2M22383	Useable

Note: If a Logic PCB with a status “not useable” is found this PCB must be replaced with a Logic PCB with a status “useable”.

**Devices affected:** All devices Radiant Heater RH600

**When:** Next repair or when a Logic PCB with a status “not useable” is found.

**Reporting:** Yes, by fax to Regulatory Affairs

**Additional Info:** N/A

**Cost:** N/A

**Ordering Info:** Logic PCB P/N 2M22383

**Distribution:** Dräger Service Personnel and Authorized Service Organizations.

If you have any questions, please contact Technical Support by phone at 1-800-543-5047 or by fax at 1-215-721-5789.

Dräger Medical, Inc.  
Technical Product Manager